

# BM Respirometry for Actual Oxygen Requirements in the aeration system of activated sludge processes



**SURCIS, S.L.**  
Emilio Serrano ([eserrano@surcis.com](mailto:eserrano@surcis.com))  
[www.surcis.com](http://www.surcis.com)


# 1. Yield coefficient ( $Y_H$ ) for heterotrophic biomass

First thing is to make a solution of 400 mg of sodium acetate in 1 litre of distilled water. For this solution, we must obtain (from the lab.) the actual COD value ( $COD_{ac} \approx 300$  mg/l)


We carry out an R assay in order to determine the consumed oxygen (CO)

CO: Consumed oxygen =  $\Delta O_2$  (mg/l)

We can make use of sample volume in between 50 and 100 ml. and pump speed to 2



Rs respirograma for CO determination



CO result

$$Y_{H,COD} = 1 - CO / COD_{ac} \quad \rightarrow \quad Y_{H,VSS} = Y_{H,CO} / 1,42 \quad \rightarrow \quad Y_{obs} = Y_{H,VSS} / (1 + 0.06 * SRT)$$

$Y_{H,COD}$ : Heterotrophic yield coefficient referred to oxygen demand (mg  $O_2$ /mg COD)

$Y_{H,COD}$ : Heterotrophic yield coefficient referred to MLVSS (mg VSS/mg COD)

$Y_{obs}$ : Observed yield (mg VSS/mg COD)


SRT: Sludge age (d)

## 2. Influent biodegradable COD ( $S_0$ )


The BM respirometer automatically calculates the total biodegradable COD (bCOD) of the influent wastewater by means a R test with influent wastewater to the biological reactor and endogenous sludge. To do this, it integrates of the exogenous respiration rate ( $R_s$ ) to calculate the consumed oxygen (CO) and apply the yield coefficient ( $Y_{H,COD}$ ):

$$bCOD = CO / (1 - Y_H)$$

Only in case the process has nitrification, you must add a dose of Allyl Thiourea (ATU) to the endogenous sludge (stirring and aerating it) 1/2 hour before the test (2 to 3 mg ATU / g VSS)



bCOD and  $R_s$  respirograms



Result

### 3. Effluent biodegradable COD ( $S_e$ )

$$S_e = \text{COD}_e - \text{iCOD}$$

$\text{COD}_e$ : Effluent COD (mg/L)

$\text{iCOD}$ : Inert COD (mg/L) =  $\text{COD}_o - S_o$

$\text{COD}_o$ : Influent COD (mg/L)

### 4. Actual Oxygen Requirement (AOR)

Here we have to distinguish between process without or with nitrification-denitrification.

#### 4.1. Actual Oxygen Requirement (AOR) calculation in a process without nitrification

##### Carbonaceous AOR

$$\text{AOR}_S = Q (S_o - S_e) / 1000 - 1.42 * P_X$$

$Q$ : Influent flow to aerobic biological reactor ( $\text{m}^3/\text{d}$ )

$S_o - S_e$ : bCOD performance in the actual process (mg/L)

$P_X$ : Sludge production ( $\text{kg VSS}/\text{d}$ ) =  $Y_{\text{obs}} * Q (S_o - S_e) / 1000$

$$\text{AOR} = \text{AOR}_S$$

$\text{AOR}$ : Total oxygen requirement ( $\text{kg O}_2/\text{d}$ )

## 4.2. Actual Oxygen Requirement (AOR) calculation in a process with nitrification-denitrification

AOR = Carbonaceous AOR ( $AOR_S$ ) + Ammonium-nitrogen AOR ( $AOR_N$ ) – Oxygen credit from denitrification ( $AOR_{DN}$ )

AOR: Total actual oxygen requirement (kg O<sub>2</sub>/d)

### Carbonaceous AOR

$$AOR_S = Q (S_o - S_e) / 1000 - 1.42 * P_X$$

$S_o - S_e$ : bCOD performance in the actual process (mg/L)

$P_X$ : Sludge production (kg VSS/d) =  $Y_{obs} * Q (S_o - S_e) / 1000$

Q: Influent flow to aerobic biological reactor (m<sup>3</sup>/d)

### Ammonium-nitrogen AOR

$$AOR_N = Q * 4.57 * S_N / 1000$$

$S_N$ : Ammonium-nitrogen concentration to nitrify (mg N-NH<sub>4</sub>/L)

### Oxygen credit from denitrification

$$AOR_{DN} = 2.86 * Q * S_N$$

We assume that all the ammonium is converted into nitrate

SURCIS, S.L.