# Practical BM Respirometry for biological nitrogen removal in a wastewater treatment



# **BM Respirometry**



#### **Three different operation modes**

While most of the respirometers on the market offer only one operation mode, the BM respirometers have three different operation modes: OUR mode, Cyclic OUR mode, and R mode. Each mode develops different respirograms for automatic parameters including D.O., Temperature, and pH (in BM-Advance) from where specific applications can be made.

In a single batch reactor, the measuring system can work as LSS and LFS batch respirometry. The system is optimized by a one-sense membrane device, that together with a dividing plate, is able to isolate the measuring chamber and avoid bubbles against the DO sensor.

#### OUR **Cyclic OUR** R This mode is making use of the The R mode corresponds to a modified The cyclic OUR mode consists of LSS respirometry type. The OUR a progressive sequence of OUR respirometry type LFS test. The mode consists of a single test to measurements, generated from measuring system can be considered as measure the OUR and/or SOUR the DO trajectory when it a completely mixed batch reactor. In fluctuates between the DO. Low this mode, we get the important parameters (by manually setting the MLVSS concentration). It and DO. High set-points that advantage to work with a small volume were set at the start of the test. also has the option the get a of samples in order to minimize the test partial SOUR for any period time for an important package of simultaneous within the respirogram. several parameters measurement.







## Main automatic parameters in BM respirometer for the different operations modes



q: Specific COD removal rate (mg COD/ mg VSS.d)

Specific U referred to MLVSS concentration.

# Nitrification



## **Optimal nitrification conditions**

|    | рН  | 7.3 to 8,3    |  |
|----|---|---------------|--|
|    | т   | > 15 to 28 °C |  |
| 1. | OD  | 1 a 3 ppm     |  |
|    | COD/TKN                                     | < 5           |  |
|    | Enough HRT for nitrification in the reactor |               |  |
|    | Without any inhibitory or toxic compound    |               |  |

#### 2. Coherent SRT & MLSS vs Temperature:



# Possible reasons from which the actual nitrifier biomass is lower than the reference value in the standard table

1. The process is not operating under one or more correct conditions of DO, pH, Temperature.

2. BOD/TKN > 5

3. BOD/TKN  $\geq$  5 + Low temperature (< 15°C) for a long time.

3. Low temperature (< 15°C) for a long period of time.,

4. Presence of toxic compound in the influent wastewater.

5. Nutrients defficiency

6. Others.

## **Actual nitrification rate**

To get the nitrification rate, we carry out a R test with ammonium chloride on equivalent concentration until reaching the maximum exogenous respiration rate. Then, Allyl Thiourea (ATU) could be added to inhibit the nitrifiers bacteria and get the endogenous sludge ready for organic matter tests where only the heterotrophic biomass is acting.



Rs respirogram after adding the NH<sub>4</sub>Cl dose

 $AUR = AUR_{max} * [DO / (K_{OA} + DO)]$ 

AUR: Actual nitrification rate (mg /L.h NH<sub>4</sub>-N)

AUR<sub>max</sub>: Nitrification rate for max. DO - 4,57: total mg O<sub>2</sub> consumed per mg NH<sub>4</sub> - N

 $Rs_{\mbox{\scriptsize N}}$ : Maximum value of exogenous respiration rate (Rs) due to nitrification

DO: Average bulk dissolved oxygen in the actual nitrification process

 $K_{OA}$ : Half saturation coefficient = 0.5 mg/l - Default value ASM2, ASM3, Henze et al 2000 - [when DO  $\ge$  3 mg/L  $\rightarrow$   $K_{OA}$  = 0]

### **Active nitrifier biomass concentration**

From the nitrification rate of an already existing nitrification process

**X**<sub>A</sub> = Y<sub>A</sub> \* 24 \* AUR \* SRT

X<sub>A</sub>: Nitrifier biomass concentration (mg/l)

 $Y_A$ : Autotrophic yield coefficient  $\approx 0.12$  - Metcalf & Eddy

SRT: Actual sludge age on which the process is operating (d)

From standard table

This table can be applied for a process that is running without any inhibition problems, under an average temperature > 20 $^{\circ}$ C, pH in between 7 and 8.5 and DO ≥ 2 ppm.



Source: Metcalf & Eddy. 1995

COD: Elimnated COD (mg/L) TKN: Eliminated TKN (mg N/L)

 $X_A = F_N * MLVSS$ 

#### Nitrogen for nitrification

Actual nitrogen being nitrified

$$\mathbf{S}_{n} = \mathsf{TKN}_{O} - \mathsf{N}_{O} - \mathsf{S}_{sy} - \mathsf{S}_{ne}$$

 $S_N$ : Actual ammonium-nitrogen concentration that is currently being nitrified (mg/l NH<sub>4</sub>-N) by the actual AUR TKN<sub>0</sub>: Influent TKN (mg/L N) N<sub>0</sub>: Organic nitrogen in the effluent  $\approx$  2 mg N/L N<sub>sy</sub>: Nitrogen utilized in the cell synthesis = 0.05 \* BOD  $S_{Ne}$ : Actual ammonium-nitrogen concentration in the effluent (mg N-NH<sub>4</sub>/L)

$$E_n = 100 * (1 - S_{ne} / S_n)$$

E<sub>n</sub>: Actual efficiency (%)

#### Required Nitrogen to be nitrified

 $\mathbf{S'_N} = \mathsf{TKN}_{\mathsf{O}} - \mathsf{N}_{\mathsf{O}} - \mathsf{S}_{\mathsf{sy}} - \mathsf{S'}_{\mathsf{ne}}$ 

 $S'_{N}$ : Ammonium-Nitrogen concentration that should be nitrified (mg NH<sub>4</sub>-N/L)  $S'_{Ne}$ : Required ammonium-nitrogen concentration in the effluent (mg N-NH<sub>4</sub>/L)

 $E'_{n} = 100 * (1 - S'_{ne} / S_{n})$ 

E'n: Required efficiency (%)

Estimation of required nitrification rate

**AUR'** =  $[1 + (E'_n - E_n) / 100] * AUR$ 

AUR': Nitrification rate necessary to get the E'<sub>n</sub> efficiency (mg NH<sub>4</sub>-N /l.h)

#### Influence of dissolved oxygen on the nitrification rate

For different DO values , the actual AUR value can be plotted.

 $AUR = AUR_{max} * DO / (K_{OA} + DO)$ 



### Minimum oxygen for any required nitrification rate

From the curve AUR vs DO, the minimum DO ( $DO_{min}$ ), on which the nitrification could operate, can be figured out for any required nitrification rate AUR'.



 $DO_{min} = K_{OA} / (AUR_{max} / AUR' - 1)$ 

### Influence of Temperature and pH on the nitrification rate

Nitrification rate depends of the corresponding exogenous respiration rate (Rs) as ammonia is being eliminated. Besides de oxygen influence, the temperature and pH have also a direct influence in the nitrification rate and therefore the knowledge of this should represent an essential data to get the actual nitrification rate.



In a BM Advance respirometer there is the option of be varying the pH during the test to analyze the different Rs and determine the break-level at which the inhibition starts

The pH influence can be analyzed on different temperature settings and it also could be analyzed the Temperature influence for a fixed pH level.

# Possible reasons for which the actual nitrification rate could be less than the required nitrification rate

- 1. The process is not operating under one or more correct conditions of DO, pH, Temperature.
- 2. The sludge age (SRT) is lower than the one on which the process should be operating.
- 3. The concentration of the active nitrifier biomass is too low.
  - This could be because of
    - . BOD/TKN > 5
    - . Conditions are out of normal range
    - . Inhibition effect
- 5. Toxicity
- 6. Others

#### Actual nitrifier biomass growing rate

 $\mu_{A}$ = 24 \* AUR \* Y<sub>A</sub> / X<sub>A</sub>

 $\mu_A$ : Nitrifier biomass growing rate (d<sup>-1</sup>)

 $Y_A$ : Yield coefficient  $\approx 0,12$  (Metcalf & Eddy)

#### Minimum sludge age for the actual nitrification

**SRT**<sub>min</sub> = 1 / ( $\mu_A - b_A$ )

SRT<sub>min</sub>: Minimum sludge age for nitrification (d)

| Temp | Death &<br>Decay Rate<br>b <sub>A</sub> (days <sup>-1</sup> ) |
|------|---|
| 10°C | 0.02  |
| 15°C | 0.03  |
| 20°C | 0.04  |
| 25°C | 0.05  |

Source: EPA – Long Island. Sound Study, NY - 2000

Sludge age condition for nitrification  $\rightarrow$  SRT  $\geq$  SRT<sub>min</sub>

# Practical operation protocol for the nitrification under energy optimization frame

In this protocol we asume the temperature and pH conditions are within the normal range and the process has the capability to control the DO level under the approach of its minimum range.



## Denitrification



## **Conditions for denitrification process**

| Conditions  |                   |   |  |  |  |
|---|-------------------|---|--|--|--|
| рН  | 6.5 a 8 (optimal) |   |  |  |  |
| BOD/TKN   | 2.5 to 5          |   |  |  |  |
| Soluble biodegradable COD/N-NO <sub>3.DN</sub>              | ≥ 2.83            | - |  |  |  |
| DO  | < 0.3 mg/L        | - |  |  |  |
| Denitrification zone with enough HRT to perform the process |                   |   |  |  |  |
| Without any inhibitor nor toxic compounds                   |                   |   |  |  |  |

### Determination of NUR by sodium acetate solution

Just as a thmb rule, the value of the estimated NUR could be determined by using the sodium acetate solution and taking the average of the COD utilization rate (U) in the R test.

Here, in order to make the corresponding COD solution, it should be taken into account that 1 mg sodium-acetate/l correspond approx. to a readily biodegradable COD ( $COD_{ac}$ ) of 0,75 mg/l.



CO Respirogram – sodium acetate solution sample

**NUR** = (U  $(1 - Y_H) / 2,86) * K'_O / (K'_O + OD_{DN})$ 

U: COD utilization rate (mg COD/l.h)

 $Y_{H}$ : Yield coefficient (mgO<sub>2</sub>/COD)  $\approx$  0,63 (defaul value)

 $K'_{o}$ : Inhibition coefficient due to oxygen in the anoxic zone = 0,2 (mg/L) - Default value –

 $DO_{DN}$ : Dissolved oxygen in the denitrification zone (mg /L  $O_2$ )

#### rbCOD required for denitrification

When the anoxic denitrification process takes place just before the aerobic treatment, a possible critical situation could arise when there is very little rbCOD in the influent and it is completely removed in that zone. Thus, there would be no rbCOD left for the aerobic zone and only the slowly biodegradable COD (rbCOD) would be in the aerobic zone. Then, under that condition and due to of a possible lack of direct food, it may result in low biomass heterotrophic growth, weak flocculation and poor process performance.

For this reason, in many cases, it is important to exceed the minimum rbCOD required for a given nitrate concentration.



#### rbCOD required for denitrification

$$rbCOD_{DN} > 2,86 * S_{NO3} / (1 - Y_{H})$$

 $rbCOD_{DN}$  Minimum readily biodegradable COD required for denitrification (mg/L) Y<sub>H</sub>: Heterotrophic yield coefficient (mg O<sub>2</sub>/mg COD)  $\approx$  0,63 (usual default value for the S<sub>NO3</sub>: Nitrate to denitrify (mg NO<sub>3</sub>-N/L)

If the required rbCOD is not met, it may be necessary to resort to the use of an external source of of readily biodegradable COD (normally methanol)

#### Calculation of methanol loading for a determined denitrification efficiency

#### PROCEDURE

- 1. Calculation of the  $S_{NO3}$  and rbCOD present in the influent
  - Calculation of the nitrate concentration to denitrify (S<sub>NO3</sub>)
  - Calculation of the rbCOD in the influent by mean a R test with filtered wastewater sample.
- 2. Calculation of the methanol loading for maximum denitrification
  - Methanol loading (m<sup>3</sup>/d) = 2.86 \* Nitrate loading rbCOD loading Nitrate loading (kg S<sub>NO3</sub>/d) = Q (m<sup>3</sup>/d) \* S<sub>NO3</sub> (kg NO<sub>3</sub>-N)/m<sup>3</sup> rbCOD loading (kg COD/d) = = Q (m<sup>3</sup>/d) \* rbCOD (kg COD/m<sup>3</sup>)
  - Methanol loading (L/d) = 1000 \* Methanol loading (m<sup>3</sup>/d)
    1 liter methanol = 1.2 Kg rbCOD
  - Methanol loading (kg/d) = Methanol loading (L/d) / 1.2
- 3. Calculation of the methanol loading for a determined denitrification efficiency (E)
  - E = 1 S<sub>NO3e</sub> / S<sub>NO3</sub>
    E: denitrification efficiency
    S<sub>NO3e</sub> : Nitrate expected in the effluent (mg NO<sub>3</sub>-N/L)
  - Methanol loading' = E x Methanol loading

Methanol loading': Methanol loading for E efficiency (liter/day)

In any case, the metahnol must be progressively acclimated to the sludge before applying the calculated methanol loading

#### Specific denitrification rate

**SDNR** = 24 \* NUR /  $X_v$ 

SDNR: Specific denitrification rate (mg N-NO<sub>3</sub> / mg VSS.d)

#### Guide table

#### Estimated Specific Denitrification Rates

| Temp<br>° C | Estimated<br>SDNR | Temp<br>° C | Estimated<br>SDNR |
|-------------|-------------------|-------------|-------------------|
| 10          | 0.035             | 18          | 0.076             |
| 12          | 0.042             | 20          | 0.091             |
| 14          | 0.052             | 22          | 0.110             |
| 16          | 0.063             | 24          | 0.132             |

Source: Long Island Sound Training – Nitrogen Removal - 2003 (EPA)

# Possible reasons for which the actual specific denitrification rate could be less than the reference table value

- 1. The process is not operating under one or more correct conditions range.
- 2. The concentration of the readily biodegradable COD is too low (soluble organic cabonaceous matter)
- 3. Anoxic zone is not gathering the anoxic condition (oxygen < 0.3 ppm)
- 4. The hidraulc retention time in the denitrification zone is too short (it has not enough volume)
- 5. Presence of inhibitor or toxic compounds in the wastewater.
- 6. Others

#### SURCIS, S.L.

#### Teléfono: 932 194 595 / 652 803 255

E-mail: <u>surcis@surcis.com</u> / <u>eserrano@surcis.com</u>

Internet: www.surcis.com